

GENESIS

Anaesthesia Workstation

SERVICE MANUAL

Training document



Document No. 9180021

Rev. 3.x.y.2 (2023-03)

Valid for software versions 2.x.y – 3.x.y

The Field Service Certification designates exclusively the **Recipient** person to conduct the procedures addressed in this documentation.

Copyright 2023 by HERSILL SL.

This document is personal and non-transferable.

Any disclosure of this document to third parties must be expressly authorized in writing by Hersill. The recipient is responsible for any unauthorized dissemination of the document.

All rights reserved. No part of this publication may be reproduced, stored in, or introduced into a retrieval system, or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise) without prior written permission of Hersill.

This manual may be revised or replaced by Hersill at any time and without notice. Ensure that you have the most current applicable version of this manual; if in doubt, contact Hersill.

While the information set forth is believed to be accurate, it is not a substitute for the exercise of professional judgment.

Nothing in this manual shall limit or restrict in any way Hersill's right to revise or otherwise change or modify the equipment (including its software) without notice. In the absence of an express, written agreement to the contrary, Hersill has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including software) described herein.

Edited and printed in Spain.

Contents

1. Symbols	7
2. Abbreviations	9
3. User and Patient Safety	10
4. Installation	11
4.1. Configuration of non-standard passive AGSS set-up	11
5. Cleaning, disinfection and sterilization	14
6. Maintenance	14
7. Calibration	15
7.1. Introduction to calibration conditions	17
7.2. Supply Pressures	18
7.3. Airway Flow	20
7.4. Airway Pressure	23
7.5. PEEP Sensor	25
7.6. Air Flow Fresh Gas	26
7.7. O2 Flow Fresh Gas	27
7.8. N2O Flow Fresh Gas	28
7.9. O2 Safety Backup Flow Fresh Gas	29
7.10. Fresh Gas Pressure	30
7.11. AGSS flow	33
7.12. Battery charge	34
8. Access levels	36
9. Execution modes	38
10. Machine setup	39
11. SelfTest	40
12. Service control	41
13. Service panel	48
13.1. Screen selection bar	48
13.2. Calibrations panel	48
13.3. Databases panel	49
13.4. Commands panel	49
13.5. Keyboard	49
13.6. Ventilation screen	50
13.7. PIDs screen	52
13.8. Inputs screen	55
13.9. Outputs screen	58
14. Diagrams	60

14.1.	General electrical diagram	60
14.2.	Control module - electrical diagram	61
14.3.	Gas manifold – electrical diagram.....	62
14.4.	General pneumatic diagram.....	63
14.5.	Gas manifold - pneumatic diagram	64
14.6.	Breathing system - pneumatic diagram	65
15.	Troubleshooting trees	66
15.1.	Air / O2 / N2O supply pressure failure	66
15.2.	Internal flow failure	67
15.3.	Internal pressure failure	68
15.4.	VC pressure failure	69
15.5.	Airway flow failure	70
15.6.	Airway pressure failure.....	71
15.7.	Air / O2 / N2O fresh gas flow failure.....	72
15.8.	Safety O2 flow-meter failure.....	73
15.9.	Fresh gas pressure failure	74
15.10.	O2 sensor failure	75
15.11.	VC system failure	76
15.12.	Air / O2 / N2O fresh gas flow failure.....	77
15.13.	Autozero failure	78
15.14.	Ventilation System failure.....	80
15.15.	Leakage and Volume failure	88
15.16.	Driving gas switch O2-Air failure	89
15.17.	Auto-Manual switch failure	90
15.18.	Disk failure.....	91
15.19.	Battery failure	92
15.20.	Multigas analyzer failure	93
15.21.	AGSS system failure	94
15.22.	Anesthetic agent vaporizing failure	95
15.23.	CABS leakage	96
15.24.	Flow graph failure.....	97
15.25.	Flush O2+ failure	98
15.26.	Fresh-gas switch (back-front) reading failure.....	99
15.27.	Gas Manifold leakage	100
15.28.	Inverted graph failure	101
15.29.	Peripherals failure	102
15.30.	Switch-on failure.....	103
15.31.	Tidal volume I vs E failure	105
15.32.	Touch screen failure.....	107
16.	Service report	108


1. Symbols



A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.



A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the equipment or other property.



A NOTE provides additional information intended to avoid inconvenience during operation.



Consult the attached documents



Do not use lubricants neither greases in the equipment



Equipotential bonding




Defibrillation-proof type BF applied part



Operating instructions



Refer to instruction manual



Gas sample (input to the side-stream multigas analyzer)

2. Abbreviations

Abbreviation	Meaning
AGSS	Anesthetic Gas Scavenging System
APL valve	Adjustable Pressure Limiting valve
AP-VG	Adaptive Pressure -Volume Guaranteed
BPM	Breaths Per Minute
C	Compliance
C.A.B.S.	Compact Anaesthesia Breathing System
EFG	External Fresh Gas
EXP	Expiration
F	Respiratory Frequency
F/V	Flow/ Volume
FiO₂	Fraction of Inspired Oxygen
Flow-P sensor	Flow-Pressure Sensor
HLD	High Level Disinfection
HF Surgical	High Frequency electro-Surgical equipment
I and E valves	Inspiration and Expiration valves
I : E	Relation Inspiratory:Expiratory time
INSP	Inspiration
M	Manual
MAC	Minimum Alveolar Concentration
MSMA	Main Stream Multi-Gas Analyzer
P_{INSP}	Inspiratory Pressure
P_{MAX}	Maximum allowable pressure
P_{PEAK}	Peak Pressure
P_{PLAT}	Plateau Pressure
P_{AW}	Pressure Air Way
PCV	Pressure Controlled Ventilation
PEEP	Positive End Expiratory Pressure
PS	Pressure Support
R	Resistance
S	Stand-by
Sp	Spontaneous Breath
SIMV	Synchronized Intermittent Mandatory Ventilation
SSMA	Side Stream Multi- Gas Analyzer
TFT	Screen type Thin -Film Transistor
Trigger	Patient threshold to start inspiration
UPS	Uninterruptible Power Supply
VCV	Volume Controlled Ventilation
V_M	Volume Minute
V_T	Tidal Volume
V_{Te}	Tidal Volume Expired
V_{Ti}	Tidal Volume Inspired

3. User and Patient Safety



WARNING

Read carefully and pursue the User and Patient safety information contents in the Instruction for Use (IFU) documentation.

The safety information below is complementary to the IFU.



WARNING

Not for use in areas of explosion hazard:

GENESIS is neither approved for use in areas where combustible or explosive gas mixtures are likely to occur.

- Explosive anesthetics, such as ether or cyclopropane, must not be used due to the risk of fire.
- Do not introduce medicaments or substances based on flammable solvents (ex. alcohol) in the patient system. Danger of fire.
- If flammable substances are used for disinfection purposes, it is necessary to guarantee the correct evaporation.
- Do not use oil or grease.
- Do not smoke.



WARNING

Safe connection with other electrical equipment:

Electrical connections to equipment which is not listed in these Instructions for use should only be made following consultations with the respective manufacturers or an expert (i.e. should confirm the system as being in compliance with Clause 16 of IEC 60601-1:2005).



CAUTION

Strictly follow the Instructions for Use and the Service Manual:

Any use of the apparatus requires full understanding and a strict observation of these instructions.

The apparatus is only to be used for purposes specified here.



CAUTION

Maintenance:

The unit must be inspected and serviced regularly by trained service personnel according to the Maintenance Intervals

Repair and general overhaul of the apparatus may only be carried out by trained service personnel.

Only original Hersill spare parts may be used for maintenance.



CAUTION

The two side handles are intended to move the GENESIS on the floor thanks to its wheels; GENESIS is a 105 kg weight that cannot be raised from its handles. Trying to raise the equipment from these handles can damage the system.

4. Installation



WARNING

Read carefully and pursue the Installation information contents in the Instruction for Use (IFU) documentation.

4.1. Configuration of non-standard passive AGSS set-up



CAUTION

According to the standard ISO 8835-3, the AGSS system of the surgical room, required in to operate an anesthesia workstation in a safe way, must provide a suction flow between 35 and 70 L/min.

GENESIS integrates in the touch-screen an AGSS performance monitor that shows the level of the suction flow according to this standard.



CAUTION

The standard ISO 8835-3, defines AGSS system that the operation room must have. Nevertheless, there are surgical rooms that do not have an AGSS installation, and GENESIS can be configured to be operated in such non-regular installations.

The standard (correct) AGSS systems are called “Active AGSS”, and the non-regular (with risk of contamination) lack of an AGSS system are called “Passive AGSS”.

GENESIS is configured in factory set-up for regular “Active AGSS”, and can be modified (hardware and software) to operate in a “Passive AGSS” environment. The required procedure must be performed only by trained technical staff.



In order to use active AGSS with GENESIS, the scavenging is established by linking the GENESIS AGSS to the AGSS in the operating room by means of a specific tube. The efficiency of the *active* AGSS is limited by the anesthetic gas scavenging device.

A passive AGSS is usually established by a wide diameter tube which links directly from GENESIS AGSS to outside from the operation room (the free air ambient). The wider and shorter the tube is, the best; but the efficiency of the passive AGSS cannot be monitored by the GENESIS sensors and a passive AGSS is usually a source of contamination for the surgical room, due to this reason it isn't accepted by the standard ISO 8835-3.


In the factory standard configuration, when GENESIS is operating in a passive AGSS environment, the AGSS performance monitor is continuously alarming in the touch screen. It is able to modify the software configuration set-up to disable the AGSS monitoring and alarm, from the service screen access level, and at lower level, it is possible to simply disable the alarm from the user alarm configuration.

To configure GENESIS for a passive AGSS set-up, it must be followed the next steps:

- Take out the C.A.B.S. (Compact Anaesthesia Breathing System) from the GENESIS frame:



- Cover the relief hole with the valve included in the Passive AGSS kit - P/N 4180512:



Valve



Hole to cover



Screw the valve with the hand



Screw to the end without force

- Place the C.A.B.S. (Compact Anaesthesia Breathing System) on the GENESIS frame.

- Fix the sticker from the Passive AGSS kit - P/N 4180512 near by the text “AGSS”:



Now, GENESIS is ready to work in a Passive AGSS installation.

5. Cleaning, disinfection and sterilization



NOTE

Read carefully and pursue the Cleaning, disinfection and sterilization procedures contents in the Instruction for Use (IFU) documentation.



NOTE

Cleaning, disinfection and sterilization procedures don't need technical expertise to be performed, and can be performed by the technical staff from the hospital surgical area.

6. Maintenance

It is important to consider the difference between 'Maintenance' (as a planned program for preventive maintenance) and 'Service' (as a fault resolution procedure, not planned).

- Maintenance (in terms of planned maintenance or preventive maintenance): Genesis is designed to make able all the planned maintenance to be performed by the technical staff from the hospital, pursuant the procedures and plan dealt in the Instructions for Use, without the need of any 'in factory' technical training.
- Service (in term of fault resolution procedures): It is dealt in this Service Manual.

A ventilation verification test of the GENESIS Anesthesia Workstation needs to be carried out after every modification due to maintenance or service tasks.



NOTE

Additionally to the replacement of the Nomoline[®] gas sampling line when needed, ISA Sidestream multi-gas analyzer measurements must be checked every year.

Please, address to Masimo[®] Service Manual in order to carry out performance verification test of ISA Sidestream multi-gas analyzer.



NOTE

For each service intervention, it must be fulfilled a new report from chapter 16.

A copy of all the service reports must be sent to the manufacturer if a defect or a failure has to be covered by the product warranty.

7. Calibration

The software interface of GENESIS allows to calibrate all the sensors in the workstation.



NOTE

Every GENESIS is in factory calibrated, this way, during the unit installation it shouldn't be required to re-calibrate the sensors. Therefore, the calibration procedure will be only recommended when a non-calibration condition is suspected or confirmed.

Tools required for GENESIS calibration:

- Ventilator analyzer, able to measure in STPD conditions:
 - With an accuracy in flow measurements of, at least, 0.05 L/min.
 - With a range in pressure measurements of, at least, 300cmH₂O.

Best option: Fluke VT650



Alternative option: TSI certifier FA Plus



- Adapter 22F-22F
- Adapter 22M-22M
- A luer connector (3180182) + silicone tube 2x5 (3180183)




To access the GENESIS calibration screen, it is required to enter at Technical or Factory access level (see chapter 8).


From the Standby screen, press the right arrow



And press the **Calibration** option



Now you are in the sensors calibration menu screen



Proceed to calibrate the sensors that are needed to solve the trouble, pursuant the next chapters.

7.1. Introduction to calibration conditions



According to the ISO 80601-2-13 (Medical electrical equipment -- Part 2-13: Particular requirements for basic safety and essential performance of an anaesthetic workstation), 201.5.101.2 (Specifications of the flow and the gas leakage): The specifications of the flow, volume and gas leakage are expressed in STPD conditions (Standard Temperature and Pressure Dry), with the exception of the associated to the Anaesthetic breathing system, that are expressed in BTPS conditions (Body Temperature and Pressure Saturated). ⁽¹⁾

Note 1 - STPD is 101,3 kPa at operation temperature of 20°C.

Note 2 - BTPS is the local atmospheric pressure and a relative humidity of 100% at an operation temperature of 37 °C.

All the test are corrected to express in STPD or BTPS.

(1) Transcription from the ISO 80601-2-13 (Medical electrical equipment -- Particular requirements for basic safety and essential performance of an anaesthetic workstation), 201.5.101.2 (Specifications of the flow and the gas leakage):

The quantities of gas are often expressed as the volume occupied by the gas under standard conditions. an atmosphere (101.3 kPa) is generally used as normal pressure. However, several normal temperatures are used. While 0°C normal temperature is used as in physics, engineering often used either 20 °C or 21.2 °C (70 °F). In the ventilation, the gas in the lungs has an identical temperature to body temperature (approximately 37 °C) regardless gas temperature administered for anaesthesia ventilator. The volume of a given amount of gas increased by approximately 13.5%, from 0 °C to 37 °C or 5.8%, from 20 °C to 37 °C.

The gas delivery systems that supply pressurized gas to medical equipment including anaesthesia ventilators, follow technical conventions and specify the amounts of gas and flow rates in STPD conditions. This practice is followed in this international standard for all requirements related to the gas inlet.

However, anaesthesia ventilators that meet this international standard is probable that inflate the patient's lungs to values of a local atmospheric pressure range of 70 kPa to 110 kPa. Besides the gas in the lungs is always moisture saturated whatever the humidity of the administered gas from anaesthesia ventilator. With a normal temperature of 0 °C, 1L of gas referred to STPD conditions can expand the lungs 1.8L at a pressure of 70 kPa. So that values are comparable between different anaesthesia ventilators, it is essential that information for all anaesthesia ventilators are referred to the same normal conditions. Since the volume of gas and not the number of molecules which expands the lungs, BTPS conditions constitute the appropriate set of reference conditions to use.



In anaesthesia ventilators are used a variety of flow transducers. While a hot wire anemometer measures the mass flow rate of gas regardless of pressure, a pneumotachograph measures the flow of gas at the present pressure. Therefore, the necessary corrections depend on the flow transducer. When a pressure correction is required, it can be estimated appropriately. Necessary corrections are also dependent on the flow transducer placement in the anaesthetic respiratory system. The gas humidity may be zero when the flow transducer measures the inspired flow within the anaesthesia ventilator. However, when the flow transducer is placed in the Y piece, the relative humidity can be any value up to 100%. When a moisture and heat exchanger is used for humidification, the transducer output depends on its placement (distal or proximal). For a turbine type anaesthesia ventilator that is using ambient air, the air moisture drawn inside may be unknown to the anaesthesia ventilator. All these effects together inevitably introduce some errors in the conversion of the flow signal measured at BTPS conditions. However, these errors are only a few percent.

7.2. Supply Pressures



Calibrate the supply pressures sensors in STPD conditions (see chapter 7.1).

Proceed according to the steps and instructions described in the screen:



1- Disconnect the three main hoses (Air, Oxygen and N₂O)





2- Disconnect the three main hoses (Air, Oxygen and N₂O)





3- Connect the air hose and introduce the pressure



4- Connect the N₂O hose and introduce the pressure



5- Connect the oxygen hose and introduce the pressure



7.3. Airway Flow

Calibrate airway flow in STPD conditions (see chapter 7.1).


10:29
27/01/16

Lamp

Adult Flow-P and Driving sensors


Help ?


↑

↑

↓

↓

Cal. P.	A.W. Ad.	Driving	-
-90.0	2490 mv	-1926 mv	
-70.0	2920 mv	-955 mv	
-50.0	3274 mv	-179 mv	
-30.0	3527 mv	391 mv	
-10.0	3674 mv	713 mv	
-5.0	3692 mv	756 mv	
0.0	3699 mv	757 mv	
5.0	3710 mv	757 mv	
10.0	3729 mv	800 mv	
30.0	3872 mv	1122 mv	
50.0	4127 mv	1693 mv	
70.0	4470 mv	2469 mv	
90.0	4898 mv	3440 mv	



CALIBRATING

- 1.- Use a calibrated STP flow equipment (reference) with a range between -100 and 100 L/min (with minimal accuracy of 0,5 L/min).
- 2.- Connect the Air supply inlet to a supply of Air pressurized between 2,7 and 6,9 bar.
- 3.- Set the zero reference on the calibrated equipment.
- 4.- Connect the flow sensor in inverted direction with the PATIENT side connected to DRIVING OUTPUT of anesthesia station and the other side to the pattern.
- 5.- Check that the calibrate equipment indicates a measure between -0,5 and 0,5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 6.- Press Next.

O₂

N₂O

Air



Next

Cancel


10:30
27/01/16

Lamp

Pediatric Flow-P sensors


Help ?


↑

↑

↓

↓

Cal. P.	A.W. Ped.	-	-
-45.0	841 mv		
-30.0	2393 mv		
-20.0	3044 mv		
-10.0	3490 mv		
-5.0	3641 mv		
-2.5	3680 mv		
0.0	3700 mv		
2.5	3715 mv		
5.0	3761 mv		
10.0	3924 mv		
20.0	4485 mv		
30.0	5439 mv		
45.0	7109 mv		



CALIBRATING

- 1.- Use a calibrated STP flow equipment (reference) with a range between -100 and 100 L/min (with minimal accuracy of 0,5 L/min).
- 2.- Connect the Air supply inlet to a supply of Air pressurized between 2,7 and 6,9 bar.
- 3.- Set the zero reference on the calibrated equipment.
- 4.- Connect the flow sensor in inverted direction with the PATIENT side connected to DRIVING OUTPUT of anesthesia station and the other side to the pattern.
- 5.- Check that the calibrate equipment indicates a measure between -0,5 and 0,5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 6.- Press Next.

O₂

N₂O

Air



Next


Cancel


HERSILL training document Recipient: Rui Carlos Dias Cordeiro Da Mata - Hemoangola

20

Proceed according to the steps described in the screen:

- 1- Connect the calibration kit as it's shown on these pictures. Make sure the arrow of the flow sensor points the anaesthesia machine instead of the flow analyzer until you reach the 0 L/min step. Once the 0L/min step is reached, change the direction of the flow-P sensor.





2- Repeat the procedure for the pediatric Flow-P sensor (Pediatric-flow-P calibration)


7.4. Airway Pressure

Calibrate airway pressure in STPD conditions (see chapter 7.1).

[illegible]

Proceed according to the steps described in the screen:

- 1- Connect the calibration kit as it's shown on these pictures. Make sure you use the adult flow-p sensor. Use the test lung to pressurize the line. Use a calibrated pressure equipment with a 0 to 110cmH₂O range.



7.5. PEEP Sensor

Calibrate PEEP pressure in STPD conditions (see chapter 7.1).

09:28
20/10/14

Peep sensor

Logs
Options
Silence
Help ?


↑

↑

↓

↓

Cal. Point	Peep	-	-
0.0	2853 mv		
30.0	2916 mv		
60.0	2978 mv		
120.0	3099 mv		
180.0	3228 mv		
240.0	3348 mv		
300.0	3479 mv		




CALIBRATING

- 1.- Use a calibrated pressure equipment with a range between -1 and 500 cmH₂O (with accuracy of 1 cm H₂O).
- 2.- Connect the Air supply input to a supply of Air or O₂ pressurized between 2.7 and 6.9 bar.
- 3.- Set the zero reference on the calibrated equipment.
- 4.- Connect the VC. Outlet of the gas manifold to the pressure inlet in the calibrated equipment by a 'T' piece and to the VC Inlet of the C.A.B.S.
- 5.- Check that the calibrate equipment indicates a measure between -1 and 1 cmH₂O. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 6.- Press OK.

Ok
Exit

Proceed according to the steps described in the screen:

Be sure to use a calibrated equipment that can measure over 300cmH₂O.



7.6. Air Flow Fresh Gas

Calibrate Air flow in STPD conditions (see chapter 7.1).

09:29
20/10/14

Air flow sensor (Fresh Gas)

Logs
Options
Silence
Help ?


↑

↑

↓

↓

Cal. Point	F.G. Air	-	-
0.0	3330 mv		
0.1	3670 mv		
0.2	3981 mv		
0.4	4482 mv		
0.7	4891 mv		
1.0	5287 mv		
1.5	5688 mv		
2.0	5988 mv		
3.0	6430 mv		
4.0	6849 mv		
8.0	7600 mv		
12.0	8222 mv		
15.0	8530 mv		




CALIBRATING

- 1.- Set the frontal switch to E.F.G. position.
- 2.- Use a calibrated STP Air flow equipment with a range between 0 and 20 L/min (with accuracy of 0.5 L/min).
- 3.- Connect the Air supply input to a supply of Air pressurized between 2.7 and 6.9 bar.
- 4.- Set the zero reference on the calibrated equipment.
- 5.- Connect the FRESH GAS Outlet of the anesthesia station to the flow input of the reference equipment.
- 6.- Check that the calibrate equipment indicates a measure between -0.5 and 0.5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 7.- Press OK.

Ok

Exit

Proceed according to the steps described in the screen:



7.7. O2 Flow Fresh Gas

Calibrate O2 flow in STPD conditions (see chapter 7.1).

09:29
20/10/14

O2 flow sensor (Fresh Gas)

CABS

Logs


Options

Silence

Help ?

CALIBRATING

Cal. Point	F.G. O2	-	-
0.0	3560 mv		
0.1	3902 mv		
0.2	4389 mv		
0.4	5131 mv		
0.7	5358 mv		
1.0	5631 mv		
1.5	6242 mv		
2.0	6649 mv		
3.0	7178 mv		
4.0	7651 mv		
8.0	8559 mv		
12.0	9062 mv		
15.0	9240 mv		



- 1.- Set the frontal switch to E.F.G. position.
- 2.- Use a calibrated STP O2 flow equipment with a range between 0 and 20 L/min (with accuracy of 0.5 L/min).
- 3.- Connect the O2 supply input to a supply of O2 pressurized between 2.7 and 6.9 bar.
- 4.- Set the zero reference on the calibrated equipment.
- 5.- Connect the FRESH GAS Outlet of the anesthesia station to the flow input of the reference equipment.
- 6.- Check that the calibrate equipment indicates a measure between -0.5 and 0.5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 7.- Press OK.

OkExit

Keep the installation prepared previously for Air, and repeat the process with O2 instead of Air.

7.8. N2O Flow Fresh Gas

Calibrate N2O flow in STPD conditions (see chapter 7.1).

09:29
20/10/14

N2O flow sensor (Fresh Gas)

CASS

Logs


Options

Silence

Help ?

CALIBRATING

Cal. Point	F.G. N2O	-	-
0.0	3225 mv		
0.1	3834 mv		
0.2	4451 mv		
0.4	4753 mv		
0.7	5299 mv		
1.0	5803 mv		
1.5	6299 mv		
2.0	6802 mv		
3.0	7258 mv		
4.0	7559 mv		
8.0	8460 mv		
12.0	8940 mv		
15.0	9182 mv		



- 1.- Set the frontal switch to E.F.G. position.
- 2.- Use a calibrated STP N2O flow equipment with a range between 0 and 20 L/min (with accuracy of 0.5 L/min).
- 3.- Connect the N2O supply input to a supply of N2O pressurized between 2.7 and 6.9 bar.
- 4.- Set the zero reference on the calibrated equipment.
- 5.- Connect the FRESH GAS Outlet of the anesthesia station to the flow input of the reference equipment.
- 6.- Check that the calibrate equipment indicates a measure between -0.5 and 0.5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 7.- Press OK.

Ok

Exit

Keep the installation prepared previously for Air, and repeat the process with N2O instead of Air.

7.9. O2 Safety Backup Flow Fresh Gas

Calibrate O2 backup flow in STPD conditions (see chapter 7.1).

09:29
20/10/14

O2 backup flow sensor (Fresh Gas)

Logs

Options

Silence

Help
?

↑

↑

↓

↓

Cal. Point	O2B F.G.	-	-
0.0	2931 mv		
0.2	3413 mv		
0.4	3636 mv		
0.7	4137 mv		
1.0	4702 mv		
1.5	5196 mv		
2.0	5693 mv		
3.0	6024 mv		
4.0	6418 mv		
6.0	6807 mv		
8.0	7240 mv		
12.0	7453 mv		
15.0	7662 mv		

CALIBRATING

- 1.- Set the frontal switch to E.F.G. position.
- 2.- Use a calibrated STP O2 flow equipment with a range between 0 and 20 L/min (with accuracy of 0.5 L/min).
- 3.- Connect the O2 supply input to a supply of O2 pressurized between 2.7 and 6.9 bar.
- 4.- Set the zero reference on the calibrated equipment.
- 5.- Connect the FRESH GAS Outlet of the anesthesia station to the flow input of the reference equipment.
- 6.- Check that the calibrate equipment indicates a measure between -0.5 and 0.5 L/min. Otherwise, contact with Technical Service (there would be an internal leak in the ventilator valve).
- 7.- Press OK.

Ok

Exit

The procedure is the same as per of O₂ Flow Fresh Gas (chapter 7.7), but to control de flow it is required to turn the manual safety O₂ knob.





7.10. Fresh Gas Pressure

Calibrate fresh gas pressure in STPD conditions (see chapter 7.1).

[illegible]

Proceed according to the steps described below.





First step: Safety O₂= 0L/min; APL valve= 0 cm H₂O.

Second step: Safety O₂= 1L/min; APL valve= 50 cm H₂O.

Introduce the value of the pressure in the keyboard.

11:15
21/10/14

Fresh gas pressure sensor

Help ?

CALIBRATING

Cal. Point	F.G.P.	-	-
0.0	29 mv		
30.0	1094 mv		

1.- Open the O2 Backup flowmeter to 1 L/min and set 50 cmH₂O in the APL.
2.- Enter the value in cmH₂O using the keyboard.
3.- Press Ok.

20.0

7 8 9
4 5 6
1 2 3
0 . ←

70.0
20.0

Accept

Test Lung

Manual Bag


Pressure Reference Instrument

Ok

Exit

7.11. AGSS flow

Proceed according to the steps described in the screen (see below).




7.12. Battery charge

See Service Panel chapter first.

- 1.- Turn off the machine, disconnect the external power and battery.
- 2.- Connect the battery cables to an external power source.
- 3.- Set the source voltage at 12 V and the current at 6 A.



- 4.- Turn on the machine again, open the service panel, press "Start Hw" and "Inputs" tab.
- 5.- In "Outputs" activate outputs 1,3 and 12.
- 6.- In "Ventilation" set a filtering with a coefficient of 0.99.



- 6.- Select the input "BAT (112)". The default calibration will appear in the right column.
- 7.- The first calibration point is 0 % and 4500 mV. That is the input voltage, of the battery sensor, at which the machine is turned off for safety. It must be modified provisionally so that it does not happen. Change the 4500 mV to 2500 mV.
- 8.- In the box on the right "Calibrations", click in this order:
 - "Make factory calibration file"
Wait for the message that you're done to appear.
 - "Recalculate calibrations in Service panel"
Wait for the message that you're done to appear.
 - "Update databases and recalculate calibration in ventilation"
Wait for the message that you're done to appear.

9.- Write the value which is read with 12 V (measured with the polimeter at the battery terminals on the machine) in "BAT 12" (100% charge).

10.- Lower the voltage of the source to 11.3 V and also point the value (0% load).



11.- Raise the source back to 12 V.

12.- Change the calibration manually in the right column. It will be necessary to enter the values pointed for 0 % and 100 % of load.

13.- In the box on the right "Calibrations", click in this order:

- "Make factory calibration file"
Wait for the message that you're done to appear.
- "Recalculate calibrations in Service panel"
Wait for the message that you're done to appear.
- "Update databases and recalculate calibration in ventilation"
Wait for the message that you're done to appear.

14.- Exit the service panel and check the calibration by seeing the percentage of battery in the power frame. With 12 V you should dial 100% and go down to 5% which should be ~11.35 V (Do not go down to 11.3 V because it will turn off).

Keep in mind that the displayed load value is very filtered, so it will drop very slowly when changing the voltage of the source.

8. Access levels

To change the access level, touch the control **[Management]** from the Standby screen.



NOTE

After Starting Up GENESIS, it is always in access level User, this way when the access level is changed, simply Shooting Down GENESIS the unit will change to User access level.

There are 4 access levels: User, Installer, Technical and Factory; they require a password. Each one allows deeper access to control the software features (see the table below):


Access level	User	Installer	Technical	Factory
password	1111	4532	2676	
Execution mode Standard	0	0	0	0
Execution mode Demo		0	0	0
Machine setup		0	0	0
Full machine setup			0	0
Execution mode Service + Service Control			0	0
See values in Self-Test			0	0
See values in Alarms			0	0
Modify the alarm preset values by default			0	0
Reset alarm preset values by default to factory values			0	0
Calibration			0	0
Technical Service Panel			0	0
Factory Service Panel				0

9. Execution modes

There are three execution modes: Standard, Demo and Service.

- Standard mode is the regular.
- Demo mode is intended basically for installing and training purposes, GENESIS substitute the patient circuit by a software simulator, and all the settings and graphics configuration are available.
- Service mode adds to the standard mode some debugging curves, in the graphics and in the controls, useful to analyze the behavior of the workstation.

		Execution mode		
		Standard	Demo	Service
Access level	User	0		
	Installer	0	0	
	Technical	0	0	0
	Factory	0	0	0



10. Machine setup

The machine setup is accessible from the Technical and Factory access levels.



NOTE

To maintain correctly the machine setup is necessary for the correct operation of the workstation. The setup is entered in factory at the end of the assembly procedure, and this setup must be updated if any module is changed after a repair procedure.

In the setup screen, it is indicated the version of each hardware module, the options enabled in the unit (advanced ventilation modes and advanced monitoring) and other special options defined by the manufacturer.

11. SelfTest

10:41
27/01/16

Basic Quick SelfTest

Help
?

The SelfTest checks the optimal operation of the ventilator and the anesthesia circuit, and helps to detect and correct faults, although the results will not be registered in any anesthesia report, leaving in the hands of the physician the decision to accept the system operation even with failed test results.

Caution!
Make sure that no patient is connected to the machine during SelfTest

Average duration of the test: 40 seconds

	Test	Status	Information	E
1	Air supply pressure		4.2 bar / 5553 mv	
2	O2 supply pressure		4.2 bar / 5782 mv	
3	N2O supply pressure		4.2 bar / 5777 mv	
4	Internal flow		0.0 L/min / 757 mv	
5	Internal pressure		0.0 cmH2O / 923 mv	
6	VC pressure		0.0 cmH2O / 2990 mv	
7	Airway flow		0.0 L/min / 3699 mv	
8	Airway pressure		0.0 cmH2O / 2047 mv	
9	Air fresh gas flow		0.0 L/min / 3594 mv	
10	O2 fresh gas flow		0.0 L/min / 3469 mv	
11	N2O fresh gas flow		0.0 L/min / 3928 mv	
12	'Safety O2' flow		0.0 L/min / 3386 mv	
13	Fresh gas pressure		1.1 cmH2O / 1741 mv	
14	O2 Concentration		0.0 % / 3999 mv	
15	AutoZero		3701.0 > 3701.0 > 3701.0 > 3701.0 > 3701.0 > 3701.0 mv	
16	Air Fresh gas supply		3591 > 3591 mv	
17	O2 Fresh gas supply		3467 > 3467 mv	
18	N2O Fresh gas supply		3928 > 3931 mv	
19	VC system		2988 > 2989 > 2989 mv	
20	Ventilation system		0 / 0 cmH2O / 0 / 0 ml	
21	Leaks and volumes		0.0 > 0.0 cmH2O	
22	Driving gas switch		0.0 > 0.0 cmH2O	
23	Auto/Manual switch		0.0 cmH2O	
24	Disk			
25	Battery			
26	Gas monitor			
27	AGS System			

Cancel



NOTE


To solve any failure in the SelfTest, apply the corresponding procedure detailed in troubleshooting chapter 15.

12. Service control

The Service control is only available from the execution mode *Service* (accessible for technical and factory access levels).

This control is useful in order to find out and solve any problem on the machine, at engineering level technical service.

Once the engineer/technician has selected the *execution mode service* the ventilation screen will change. It will go on showing the same number of tables defined on the screen graphics design (maximum 3G/3C), but in each of these graphs it can be displayed many different parameters.



- **Flow graph**

- Adult FAW
Curve of airway flow for adult Flow-P in units of flow.
- Driving Flow
Curve of driving (internal flow) in units of flow.
- Target Flow
Curve of objective flow in units of flow.
- Voltage Target
Curve of objective flow in millivolts.
- V. Target Amp. Ctr.
Curve of objective flow with amplitude control in millivolts.
- DC Input
Curve of driving control (FAW or FD) input in millivolts.
- VC Input
Curve of ventilation control (PAW) input in millivolts.
- DC Output
Curve of driving control (solenoid) output in millivolts.
- VC Output
Curve of ventilation control (solenoid) output in millivolts.

- **Volume graph**


- Volume Curve of volume in units of volume.
- Driving volume Curve of volume supply by driving flow in units of volume.
- F. Gas volume Curve of volume supply by fresh gas flow in units of vol.
- Total Adm. Vol. Curve of volume supply by driving and fresh gas in units of vol.


- **Pressure graph**

- PAW Curve of airway pressure in units of pressure.
- Driving Pressure Curve of driving (internal pressure) in units of pressure.
- Target D. Pressure Curve of objective pressure in units of pressure.
- Voltage Target Curve of objective pressure in millivolts.
- V. Target Amp. Ctr. Curve of objective pres. with amplitude control in mV.
- DC Input Curve of driving control (PAW) input in millivolts.
- VC Input Curve of ventilation control (PAW) input in millivolts.
- DC Output Curve of driving control (solenoid) output in millivolts.
- VC Output Curve of ventilation control (solenoid) output in millivolts.
- VC Press. Curve of VC pressure (under the membrane) in millivolts.

Beside these options, there is a new possibility on the chart zone that is called “service control”.

When the service control is selected the technician can displayed ten different charts on the same area.





- **PID Control + A.C. Inspiration**

- P I P Primary (Driving) – Inspiration – Proportional gain and output.
- P I I Primary (Driving) – Inspiration – Integral gain and output.
- P I D Primary (Driving) – Inspiration – Derivative gain and output.
- S I P Secondary (Driving) – Inspiration – Proportional gain and output.
- S I I Secondary (Driving) – Inspiration – Integral gain and output.
- S I D Secondary (Driving) – Inspiration – Derivative gain and output.
- IIni Initial value of amplitude control.
- IAccIni Accumulated initial (init. of cycle) value of amplitude control.
- IAcc Accumulated value of amplitude control.
- Cycle Number of cycles and type of control.

- **PID Control + A.C. PEEP**

- P E P Primary (Driving) – Expiration – Proportional gain and output.
- P E I Primary (Driving) – Expiration – Integral gain and output.
- P E D Primary (Driving) – Expiration – Derivative gain and output.
- S E P Secondary (Driving) – Expiration – Proportional gain and output.
- S E I Secondary (Driving) – Expiration – Integral gain and output.
- S E D Secondary (Driving) – Expiration – Derivative gain and output.
- IIni Initial value of amplitude control.
- IAccIni Accumulated initial (init. of cycle) value of amplitude control.
- IAcc Accumulated value of amplitude control.

- **PID Control + A.C. Support**

- P E P Primary (Driving) – Expiration – Proportional gain and output.
- P E I Primary (Driving) – Expiration – Integral gain and output.
- P E D Primary (Driving) – Expiration – Derivative gain and output.
- S E P Secondary (Driving) – Expiration – Proportional gain and output.
- S E I Secondary (Driving) – Expiration – Integral gain and output.
- S E D Secondary (Driving) – Expiration – Derivative gain and output.
- IIni Initial value of amplitude control.
- IAccIni Accumulated initial (init. of cycle) value of amplitude control.
- IAcc Accumulated value of amplitude control.

- **PID Control + F. G. N2O**

- P Proportional gain and output.
- I Integral gain and output.
- D Derivative gain and output.
- Input F.G. N2O flow (sensor) in millivolts.
- Target F.G. N2O objective in millivolts.
- Output F.G. N2O output (solenoid) in millivolts.

- **PID Control + F. G. Air**

- P Proportional gain and output.
- I Integral gain and output.
- D Derivative gain and output.
- Input F.G. Air flow (sensor) in millivolts.
- Target F.G. Air objective in millivolts.
- Output F.G. Air output (solenoid) in millivolts.

- **PID Control + F. G. O2**

- P Proportional gain and output.
- I Integral gain and output.
- D Derivative gain and output.
- Input F.G. O2 flow (sensor) in millivolts.
- Target F.G. O2 objective in millivolts.
- Output F.G. O2 output (solenoid) in millivolts.

- **Auxiliary sensors**

- F N2O Fresh gas N2O flow in millivolts and in units of flow.
- F Air Fresh gas Air flow in millivolts and in units of flow.
- F O2 Fresh gas O2 flow in millivolts and in units of flow.
- F Safety O2 Fresh gas safety O2 flow in millivolts and in units of flow.
- F AGSS AGSS extraction flow in millivolts and in units of flow.
- P FG Fresh gas pressure in millivolts and in unit of pressure.
- Output F.G. N2O output (solenoid) in millivolts.

- **Pressures**

- Total Max. Maximum pressure reached in total cycle.
- Partial Max. Maximum pressure reached in partial cycle.
- Insp. C. Max. Maximum pressure reached in mandatory inspiration.
- VC Max. Maximum pressure reached in Ventilation control.
- Total Mean Mean pressure reached in total cycle.
- Partial Mean Mean pressure reached in partial cycle.
- Insp. C. Mean Mean pressure reached in mandatory inspiration.
- Total Min Minimum pressure reached in total cycle.
- Partial Min Minimum pressure reached in partial cycle.
- Insp. C. Min Minimum pressure reached in mandatory inspiration.

- **Patient volumes**

- VTi total Inspiratory tidal volume in total cycle.
- VTi partial Inspiratory tidal volume in partial cycle.
- VTi insp cycle Inspiratory tidal volume in inspiratory cycle.
- VTe total Expiratory tidal volume in total cycle.
- VTe partial Expiratory tidal volume in partial cycle.
- VTe insp cycle Expiratory tidal volume in inspiratory cycle.
- VMe total Expiratory volume minute in total cycle.
- VMe partial Expiratory volume minute in partial cycle.

- **Managed volumes**

- VTd total Driving volume supplied in total cycle.
- VTd partial Driving volume supplied in partial cycle.
- VTd insp cycle Driving volume supplied in inspiratory cycle.
- VTfg total Fresh gas volume supplied in total cycle.
- VTfg partial Fresh gas volume supplied in partial cycle.
- VTfg insp cycle Fresh gas volume supplied in inspiratory cycle.
- VTa total Driving+Fresh gas volume supplied in total cycle.
- VTa partial Driving+Fresh gas volume supplied in partial cycle.
- VTa insp cycle Driving+Fresh gas volume supplied during insp. Cycle.

- **Times of last cycle**

- D. Cycle Pr. Duration of cycle programmed.
- D. Insp. Pr. Duration of inspiration programmed.
- D. Exp. Pr. Duration of expiration programmed.
- D. Compens. Duration of compensation (Advanced cycle in SIMV).
- D. Exp. Ad. by Supp. Duration of additional exp. Time due to support.
- D. Pause Pr. Duration of inspiration pause programmed.
- D. Ramp Pr. Duration of ramp time programmed.
- T. In. C. Par. Time of partial inspiration cycle.
- D. Insp. P. Effect. Duration of partial inspiration effective.
- D. Exp. P. Effect. Duration of partial expiration effective.

○

- **Times of penultimate cycle**

- D. Cycle Pr. Duration of cycle programmed.
- D. Cycle Effect. Duration of cycle effective.
- D. Insp. Pr. Duration of inspiration programmed.
- D. Insp. P. Effect. Duration of partial inspiration effective.
- D. Exp. Pr. Duration of expiration programmed.
- D. Exp. P. Effect. Duration of partial expiration effective.
- D. Pause Pr. Duration of inspiration pause programmed.
- D. Ramp Pr. Duration of ramp time programmed.

- **Ventilation parameters**

- P Insp. Inspiration pressure programmed.
- Previous P. Insp. Inspiration pressure programmed in previous cycle.
- P. Support Support pressure programmed.
- Previous P. Support. Support pressure programmed in previous cycle.
- PEEP PEEP programmed.
- Previous PEEP PEEP programmed in previous cycle.
- Tidal volume Tidal volume programmed.
- Previous VT Tidal volume programmed in previous cycle.
- Trigger Trigger programmed.
- Previous Trigger Trigger programmed in previous cycle.

- **Gas Analyzer 1**

- Status Plugged / unplugged.
- Sensor mode SelfTest / Demo / Measure.
- O2 Compens. O2 compensation.
- N2O Compens. N2O compensation.
- Apnea time Apnea time.
- Disab. Zero Disabled zero.
- Zero in prog. Zero calibration in progress.
- Span Error Span error.
- Span cal. in prog. Span calibration in progress.

- **Gas Analyzer 2**

- Register Cfg. 0 Register of configuration 0.
- Register Cfg. 1 Register of configuration 1.
- Hw. Rev. Hardware revision.
- Sw. Rev. Software revision.
- Comm. Rev. Communication revision.
- Disab. Zero Disabled zero.
- Serial N. Serial number.

- **System**

- | | |
|-----------------------|--------------------------|
| ○ Disk free space | Disk free space. |
| ○ Disk total space | Disk total space. |
| ○ Task virtual memory | Task virtual memory. |
| ○ Task Res. Mem. | Task Res. Memory. |
| ○ % CPU User | Use of user processor. |
| ○ % CPU System | Use of system processor. |

- **Others**

- | | |
|------------------------|--|
| ○ P. Tgt. VC Insp. | Target VC pressure for inspiratory. |
| ○ P. Tgt. VC Exp. | Target VC pressure for expiratory. |
| ○ P. Tgt. AP-VG | Target AP-VG pressure for inspiratory. |
| ○ P. Tgt. AP-VG + PEEP | Target AP-VG pressure for expiratory. |
| ○ Flow Max T. | Total maximum flow. |
| ○ Flow Max P. | Partial maximum flow. |
| ○ Flow Mean I.C. | Insp. Cycle mean flow. |


13. Service panel

To access the Service panel, it is required to enter at Technical or Factory access level (see chapter 8).

From the Standby screen, press the right arrow



And press the **Service Panel** option



13.1. Screen selection bar


The bottom bar let change between the different screens of the service panel; each screen is explained forward.



13.2. Calibrations panel

Placed in the upper position from the right side panels.

To access the Calibration panel it is required to enter at Factory access level.



- Load the factory calibration file: Load in the memory the data of the factory calibration.
- Create factory calibration file: Keep safe in the calibration data the current calibration.
- Recalculate the calibration in S.P.: Recalculate the calibration using the current safe data.
- Update the Database and recalculate machine calibration: Update the data base and recalculate the calibration using the current safe data.

13.3. Databases panel


Placed in the second position from the right side panels.



- Fix register database: Repair the mistakes and the inconsistencies on the database of the register.
- Empty the register of the database: Delete all data of the database.
- Fix configuration of the database: Repair the mistakes and the inconsistencies on the database of the configuration.

13.4. Commands panel


Placed in the third position from the right side panels.



- Start HW: Activate the acquisition and reading of the PCB's data.
- Exit to console: Let close graphic interface and go to the command console.
- Enable hardware: Sends manually pulses to deactivate the timeout of the watchdog.
- Initialize PID: Initialize all the parameters to the defaults of the driver PID.
- Enable Amp.C: Activate the amplitude control in the drivers of the PID.

13.5. Keyboard

Placed in the last position from the right side panels.



- The keyboard let introduce numerical values. The up arrow let increase the selected value and the down arrow decrease it. The left arrow let delete the value.

13.6. Ventilation screen

There are four different groups of ventilation parameters on the ventilation screen: Ventilatory parameters, visualization, digital filter and units.

The screenshot displays the 'Ventilation' screen interface. At the top left is a large plot area with a y-axis ranging from 0 to 10,000 and an x-axis from 0 to 700. Below the plot are four main parameter groups:

- Ventilatory parameters:** Includes 'Mode' (0- Manual), 'Fresh Gas' (1- Air), and three flow rate settings (F.G.F. Air, F.G.F. O2, F.G.F. N2O) all set to 4000 mV.
- Unidades:** Includes 'Unit P' (10 - mV), 'Unit A' (16 - mV), 'Unit Flow' (24 - mV), 'Unit %' (1 - mV), and 'Unidad G' (8 - % Vol).
- Visualization:** Includes 'Curve 1' (1 - Consigna V), 'Curve 2' (8 - EA 1 PVA (101)), 'Curve 3' (26 - SA0 100 DG (20)), 'Curve 4' (27 - SA1 20 VC (19)), 'Min Y' (0), 'Max X' (10000), 'Filled' (0 - NO), and 'Refresh' (0,10 s).
- Digital filter:** Includes 'D. F. Kind' (0 - Basic) and 'Coefficient 1' (0,1000000000).

On the right side, there are three sections of controls:

- Calibrations:** Buttons for 'Load factory cal. file', 'Create factory cal. file', 'Recalculate calibration in S. P.', and 'Update D.B. and recal. M. cal.'.
- Database:** Buttons for 'Fix register database', 'Empty register database', and 'Fix Cfg. database'.
- Commands:** Buttons for 'Start Hw', 'Enable Hw', 'Exit to console', 'Init PID', and 'Enable Amp. C.'.

At the bottom, there is a 'Keyboard' section with a numeric keypad and directional arrows. A navigation bar at the very bottom contains buttons for 'Ventilation', 'PID Ctr. Param.', 'Inputs', 'Outputs', 'MG Analyzer', and 'Close'.

- Ventilatory parameters

- Mode

- It's used to make ventilation in pressure or volume mode. The manual mode let fix all the exit signals manually.

- Fresh Gas

- Air/ O2 /N2O

- Fresh gas flow Air

- The tension of the fresh gas flow of air can be selected

- Fresh gas flow O2

- The tension of the fresh gas flow of O2 can be selected

- Fresh gas flow N2O

- The tension of the fresh gas flow of N2O can be selected

- Visualization:

- Curve 1
Can be an analogic input or and analogic output with the unit defined.
- Curve 2
Can be an analogic input or and analogic output with the unit defined.
- Curve 3
Can be an analogic input or and analogic output with the unit defined.
- Curve 4
Can be an analogic input or and analogic output with the unit defined.
- Min Y
Minimum value of the vertical axis
- Max Y
Maximum value of the vertical axis
- Filled
The curves can be filled in the graph
- Refresh
The period of time to update the graph

- Digital filter:

There are three kinds of filter: Basic, FIR and IIR.

It is a digital filter to the analogic inputs. For the filter of analogic input, the value is a constant. (100 high filtered – 0 low filtered).

- Units

- Unit P
Units for pressure measures (mv, Pa, mmHg, HPa, mbar, cm H2O,
- Unit A
Units for pressure supply measures (mv, Pa, KPa, bar)
- Unit Flow
Units for flow measures (mv, L/min, L/s)
- unit %
Units for % measures (mv / %)
- Unit G
Units for gas concentration measures (mmHg,%vol, kPa,)

13.7. PIDs screen

On the PID Crt. Param screen there are three different groups of configuration and/or motorized signals: Control driving, Insp/ Esp. Peep control and F.G Control Air.

The screenshot displays the 'PID Crt. Param' screen with the following sections:

- Control Driving:**
 - Enable ctrl.: 0 - NO
 - Input channel: 1
 - Output channel: 0
 - Offset (mV): 3000
 - G kP: 5,000
 - G kI: 0,004
 - G kD: 10,000
 - Error P: 0,00
 - Error I: 0,00
 - Error D: 0,00
- Insp./Esp. Peep control:**
 - Enable ctrl.: 0 - NO
 - Input C.: 1
 - Output C.: 1
 - Offset (mV): 3000
 - G kP: 1,000
 - G kI: 0,002
 - G kD: 1
 - Target P. Peep (mV): 5000
 - Error P: 0,00
 - Error I: 0,00
 - Error D: 0,00
- F.G. Control Air:**
 - Enable ctrl.: 0 - NO
 - Input channel: 5
 - Output channel: 2
 - Offset (mV): 3000
 - G kP: 9,000
 - G kI: 0,100
 - G kD: 0,000
 - Error P: 0,00
 - Error I: 0,00
 - Error D: 0,00
- Calibrations:**
 - Load factory cal. file
 - Create factory cal. file
 - Recalculate calibration in S. P.
 - Update D.B. and recal. M. cal.
- Database:**
 - Fix register database
 - Empty register database
 - Fix Cfg. database
- Commands:**
 - Start Hw
 - Enable Hw
 - Exit to console
 - Init PID
 - Enable Amp. C.
- Keyboard:**
 - Numeric keypad (0-9, ., -)
 - Directional arrows (Left, Up, Down)

Bottom navigation buttons: Ventilation, PID Ctr. Param. (highlighted), Inputs, Outputs, MG Analyzer, Close.

- Control Driving
 - Enable Ctrl
 - Enable and disable the PID control of the driving
 - Input channel
 - Digital input to control PID of the driving
 - Output channel
 - Digital output to control PID of the driving
 - Offset (mv)
 - Output offset to control PID of the driving
 - G kP
 - Proportional gain for the PID of the driving
 - G kI
 - Integral gain for the PID of the driving

- G kP
Derivative gain for the PID of the driving
- Error P
Proportional error for the PID control of the driving (It's not an input parameter, it's a motorized signal)
- Error I
Integral error for the PID control of the driving (It's not an input parameter, it's a motorized signal)
- Error D
Derivative error for the PID control of the driving (It's not an input parameter, it's a motorized signal)
- Insp/Exp Peep control
 - Enable ctrl.
Enable or disable the PID control of the inspiratory and expiratory PEEP
 - Input Channel
Digital input for the PID control of the inspiratory and expiratory PEEP
 - Output channel
Digital output for the PID control of the inspiratory and expiratory PEEP
 - Offset (mv)
Offset for the PID control of the inspiratory and expiratory PEEP
 - G kP
Proportional gain for the PID control of the inspiratory and expiratory PEEP
 - G kI
Integral gain for the PID control of the inspiratory and expiratory PEEP
 - G kD
Derivative gain for the PID control of the inspiratory and expiratory PEEP
 - Target P. Peep (mV)
Target pressure in the VC chamber inspiratory and expiratory
 - Error P
Proportional error for the PID control of inspiratory and expiratory PEEP (It's not an input parameter, it's a motorized signal).
 - Error I
Integral error for the PID control of inspiratory and expiratory PEEP (It's not an input parameter, it's a motorized signal).

- Error D
Derivative error for the PID control of inspiratory and expiratory PEEP (It's not an input parameter, it's a motorized signal).

- F.G. Control (Air-O₂-N₂O)


This group of parameters changes according to the gas selected on the ventilation screen.

- Enable ctrl.
Enable or disable the PID control of the fresh gas (Air-O₂-N₂O)
- Input channel
Digital input to control the PID of the fresh gas (Air-O₂-N₂O)
- Output channel
Digital output to control the PID of the fresh gas (Air-O₂-N₂O)
- Offset (mv)
Offset for the PID control of the fresh gas (Air-O₂-N₂O)
- G kP
Proportional gain for the PID control of the fresh gas (Air-O₂-N₂O)
- G kI
Integral gain for the PID control of the fresh gas (Air-O₂-N₂O)
- G kD
Derivative gain for the PID control of the fresh gas (Air-O₂-N₂O)
- Error P
Proportional error for the PID control of the fresh gas (Air-O₂-N₂O) (It's not an input parameter, it's a motorized signal).
- Error I
Integral error for the PID control of the fresh gas (Air-O₂-N₂) (It's not an input parameter, it's a motorized signal).
- Error D
Derivative error for the PID control of the fresh gas (Air-O₂-N₂) (It's not an input parameter, it's a motorized signal).

13.8. Inputs screen

There are three different configurations in the inputs screen: Digital inputs, Analog inputs and sensor Cal.20 (Value-Voltage)

A refresh is done sequentially; it can bring about a delayed effect.



- Digital inputs:
 - On(1) / Off (0) button
Signal of the front switch button.
 - Switch back panel
Reserved for future
 - Conversation finish ADC1
Signal of conversation finish of the first converter analogic-digital.
 - Conversation finish ADC2
Signal of conversation finish of the second converter analogic-digital.
 - Power supply detection
Generated signal by the power supply source to show if it is receiving supply tension.

- Fresh gas switch gas EFG.
Signal that comes from the frontal fresh gas switch and indicates the position
- Switch fresh gas CABS.
A Signal that comes from the frontal fresh gas switch and indicates if the switch is heading CABS

The rest of the digital inputs are kept for the future.

- Analogic inputs

The analogic inputs will be show on the screen according to the configured units on the ventilation screen. There is a single column on the left side (CAL), these bottoms allow you to see the calibration values and they will be shown on the right column (Sensor CAL.) There is another column on the right side that let invert the value of each analogic input (If the tension of a sensor decreases when the value grows, it's inverted)

- DGP (110)
Pressure driving gas signal (Absolute pressure)
- PAW (101)
Pressure airway signal (Absolute pressure)
- FAW (102-100)
Flow airway signal (Differential pressure)
- O2F (26)
Oxygen fresh gas flow (Mass-flow)
- DGF (109)
Driving gas flow signal (Differential pressure)
- AIRF (25)
Air fresh gas flow signal (Mass-flow)
- BAT (112)
It indicates the tension of the battery
- VC (111)
VC – Peep pressure signal (Absolute pressure)
- O2BF (18)
Oxygen backup fresh gas flow signal (Mass-flow)

- PFG (108)
Fresh gas pressure signal (Absolute pressure)
- AGSF (122)
AGS flow signal (Differential pressure)
- RES 1 (XXX)
Reserved 1 (Differential pressure)
- RES 2 (XXX)
Reserved 2 (Differential pressure)
- RES 3 (XXX)
Reserved 3 (Differential pressure)
- N2OF (27)
Nitrous oxide flow signal (Mass-flow)
- O2P (104)
Oxygen supply pressure signal (Absolute pressure)
- N2OP (105)
Nitrous oxide supply pressure signal (Absolute pressure)
- CONC (113)
Oxygen concentration signal (Galvanic sensor)
- AIRP (103)
Air supply pressure signal (Absolute pressure)
- FAW (102-160)
Airway flow signal (Differential pressure)


- Sensor Cal

It shows the calibration chart of the sensor chosen on the CAL column.

13.9. Outputs screen

There are three different groups of parameters on the output screen. Digital outputs, analog outputs and Polynomial interpolation

A refresh is done sequentially; it can bring about a delayed effect.




- Digital outputs
 - Enable Prop. Electro valve
Activation proportional inputs signal
 - Auto-zero evs (Version depending).
The auto-zero electro valves are piloted by this signal (Depending on the version).
 - DACs Enable
Enable of the digital-analogic converters signal
 - EV Dr.Gas Air (1) / O2 (0)
The electro valve which selects the driving gas is piloted by this signal
 - BAT control sensor
This signal commutes the tension of the analog input BAT to the tension of charge.


- Turn off flag
It activates the electronic of turned off (If it's on, the machine restarts)
- Battery supply forced
It uses the battery when it is activated
- Evs, DACs and lights supply
It let feed the circuits of the electro valves, digital-analog converters and LED light of the monitor
- Multiplexer address S0
First bit of the multiplexer address of the analog inputs
- Multiplexer address S1
Second bit of the multiplexer address of the analog inputs
- Multiplexer address S2
Third bit of the multiplexer address of the analog inputs
- Gas analyzer supply
It electrically feeds the gas analyzer
- Analog outputs
 - Ev. P. Driving (20)
Piloted tension of the proportional electro valve of driving gas
 - Ev. P. VC (19)
Piloted tension of the proportional electro valve of VC
 - Ev. P. GF Aire (21)
Piloted tension of the proportional electro valve of air fresh gas
 - Ev. P. GF O2 (22)
Piloted tension of the proportional electro valve of oxygen fresh gas
 - Reserved 1
Piloted tension of the proportional electro valve of reserved 1
 - Ev. P. GF N2O (24)
Piloted tension of the proportional electro valve of nitrous oxide fresh gas
 - Reserved 2
Piloted tension of the proportional electro valve of reserved 2
 - Reserved 3
Piloted tension of the proportional electro valve of reserved 3

14. Diagrams


14.1. General electrical diagram




14.2. Control module - electrical diagram




14.3. Gas manifold – electrical diagram




14.4. General pneumatic diagram



14.5. Gas manifold - pneumatic diagram



14.6. Breathing system - pneumatic diagram



15. Troubleshooting trees



WARNING

After of finishing any troubleshooting procedure, involving hardware and parts changes or repairs, it is compulsory to proceed with the corresponding Electrical Safety Test (PICM 59.1 or according to national or local regulations).




NOTE

After of finishing any troubleshooting procedure, involving hardware and parts changes or repairs, it is needed to proceed with the corresponding calibration (see chapter 7).


15.1. Air / O₂ / N₂O supply pressure failure

Autotest 1, 2 and 3




15.2. Internal flow failure

Autotest 4




15.3. Internal pressure failure

Autotest 5




15.4. VC pressure failure

Autotest 6




15.5. Airway flow failure

Autotest 7




15.6. Airway pressure failure

Autotest 8




15.7. Air / O₂ / N₂O fresh gas flow failure

Autotest 9, 10 and 11




15.8. Safety O2 flow-meter failure

Autotest 12




15.9. Fresh gas pressure failure

Autotest 13




15.10. O2 sensor failure

Autotest 14




15.11. VC system failure

Autotest 15




15.12. Air / O2 / N2O fresh gas flow failure


Autotest 16, 17 and 18



15.13. Autozero failure


Autotest 19








15.14. Ventilation System failure


Autotest 20






















15.15. Leakage and Volume failure

Autotest 21




15.16. Driving gas switch O2-Air failure

Autotest 22




15.17. Auto-Manual switch failure

Autotest 23




15.18. Disk failure

Autotest 24




15.19. Battery failure

Autotest 25




15.20. Multigas analyzer failure

Autotest 26




15.21. AGSS system failure


Autotest 27



15.22. Anesthetic agent vaporizing failure




15.23. CABS leakage



15.24. Flow graph failure

Flow graph failure

The flow graph failure is confirmed if the Flow vs Time curve presents unusual peaks, and the peaks become bigger if the pressure increases



Conclusion:
There is a leakage in the Flow-P piping

Field Service:

- Change the two filters from the Flow-P.
- Verify / change the o-ring on the luer-lock piece for the filter connection.
- Verify / change the pair silicone tubes (circular - hexagon).

Verify / change the Flow-P sensor.


- Change the Gas manifold.

or

Factory Service:

- Take out the front plate and verify the two blue tubing.
- Take out the Gas Manifold and verify A and I ports to the pneumatic multi-contact.
 - Verify the tubing in the gas manifold.
 - Look for leaks in the auto-zero valves.

15.25. Flush O2+ failure



15.26. Fresh-gas switch (back-front) reading failure

Fresh-gas switch (back-front) reading failure

The fresh-gas switch reading failure is confirmed if changing the position of the fresh gas switch, the screen does not show the change

Field Service:


- Verify / update the Version of the Machine and the Control Unit.
- Take out the plate below the Gas Manifold (the ventilator and speaker plate).
 - Verify the cables to the Control Module.
- Take out the plate to the left side in the selectatec area.
- Unplug the Pneumatic multi-contact (to screws in the Gas Manifold)
 - Take out the Gas Manifold.
- Take out the switch unit by releasing its panel screws.
 - Change the Switch.
 - Change the Control Unit.

or


Factory Service:

- Repair the Switch: Change / move the sensors.
- Verify the cables in the Control Module: from the side connector to the PCB Control.
 - Change / Repair the PCB Control.

15.27. Gas Manifold leakage



15.28. Inverted graph failure



15.29. Peripherals failure

Peripherals failure

The peripherals failure is confirmed if there is a leak or malfunction in one of these units:

Fresh Gas Switch Back / Front
Vacuum On-Off Switch
Vacuum regulation knob
Auxiliary O2 valve
Auxiliary O2 flowmeter
Safety O2 flowmeter
Relief valve
Vacuum valve
Selectatec bar
+O2 flush

Field Service:


- Take out the plate below the Gas Manifold (the ventilator and speaker plate).
 - Take out the plate to the left side in the selectatec area.
- Unplug the Pneumatic multi-contact (to screws in the Gas Manifold)
 - Take out the Gas Manifold.
- Verify the tubing from the involved unit to the pneumatic multi-contact (position, connection, folds).
 - Take out the involved unit by releasing its panel screws.
 - Change the unit.
 - Change the gas manifold.


or

Factory Service:


- Repair the unit: If there is a leak: verify / change the o-rings; If there is a malfunction: Verify / change the damaged parts.
- Repair the Gas manifold: Verify / change the internal tubing.


15.30. Switch-on failure






15.31. Tidal volume I vs E failure





15.32. Touch screen failure



16. Service report

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

System information:	
Model: GENESIS	Serial Number:
Customer information:	
Hospital Name:	City:
Country:	Address:
Service information::	
Company:	Date:
Failure report:	Spare Parts changed (P/N): Systems changed (P/N & S/N):
Calibration equipments (Model & S/N):	Electrical Safety Test equipments (Model & S/N):
Service	Hospital
Name:	Name:
Signature:	Signature:

Manufacturer : HERSILL, S.L. – Puerto de Navacerrada 3; 28935, Móstoles, Madrid, Spain
Tel.: +34 916166000 / Fax: +34 916164892 / www.hersill.com
